



## Uv Chamber



- UV-C Spectrum we are using for better result.
- UV Intensity: 9/15/30/45/60/75 watts.
- Capacity: 18, 25, 30, 60, 75, 85 litres.
- Digital Controlled panel for better accuracy.
- Double Walled Construction.
- **Inner Chamber:** Mirror finish Stainless Steel.
- **Outer Walled** : MS with powder coated

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number** : +91 171-4005686, **Mobile**: +91-9416966084, 9896055375

**Email ID**: apexambala@gmail.com, apexambala@rediffmail.com

**Website**: [www.apexambala.com](http://www.apexambala.com)