



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Vacuum Pump (Oil Free)



Diaphragm type and Hence No Rubbing of Parts and Negligible Wear & Tear, No Lubrication Required.

(Specially For Filtration Assembly, Available In Capacities 25 Ltrs/50 Ltrs)

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)