



## Water Bath (Circulating)



- Temperature is controlled by Thermostat.
- Inner Chamber Stainless Steel & Outer MS duly powder coated/S.S.
- Size of water bath: 8" x 8" x 8".
- Temperature range: Ambient +5°C to 90°C.
- Provided with variable speed regulator to control the speed of circulation.

### Optional

- Digital temperature controller with an accuracy of  $\pm 1^{\circ}\text{C}$ .
- Microprocessor based P.I.D. based temperature controller with an accuracy of  $\pm 0.5^{\circ}\text{C}$ .

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)