



Water Bath Incubator Shaker



- Temperature range from room temperature to 100°C.
- Temp. is Controlled by Thermostat / Digital Temp. Controller / PID with an accuracy of $\pm 0.5^{\circ}\text{C}$.
- Double walled inside made of stainless steel 304 Quality and outside Mild Steel sheet painted in epoxy powder coating/S.S.
- Oscillating tray is riding on ball rollers.
- Shaking speed range 40 to 140 rpm.
- To work on 220/230 volts A.C. Supply.
- SHAKING TRAY MADE OF Stainless Steel.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Models :

(A)	275x275x75 mm.	6 Ltrs.
(B)	275x275x150mm.	12 Ltrs.
(C)	405x300x150mm	18 Ltrs.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com