



Water Bath (With Holes)



- Double walled construction with inner chamber made of Stainless Steel & outer made of S.S/Mild steel duly Powder Coated.
- S.S. Top having 75 mm holes with concentric rings.
- Temperature range upto 95°C.
- Accuracy $\pm 5^{\circ}\text{C}$ for thermostatically Controlled fitted with Neon Lamp; indicator, main switch and workable on 220 Volt 50 Hz. A.C.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com



•

Working size:-	
325X250X100 mm	6 HOLE OF 75 mm DIA
425X325X100 mm	12 HOLE OF 75 mm DIA

Optional

- Digital temperature controller with an accuracy of $\pm 1^{\circ}\text{C}$.
- Microprocessor based P.I.D. based temperature controller with an accuracy of $\pm 0.5^{\circ}\text{C}$.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com